

Author(s): MELO 3D Project Team, 2012

License: Unless otherwise noted, this material is made available under the terms of the **Creative Commons Attribution – Share-Alike 3.0 License:**
<http://creativecommons.org/licenses/by-sa/2.0/>

We have reviewed this material in accordance with U.S. Copyright Law **and have tried to maximize your ability to use, share, and adapt it.** The attribution key found at the end of the module provides information about how you may share and adapt this material.

Copyright holders of content included in this material should contact **open.michigan@umich.edu** with any questions, corrections, or clarification regarding the use of content.

For more information about **how to cite** these materials visit <http://open.umich.edu/education/about/terms-of-use>.

Any **medical information** in this material is intended to inform and educate and is not a tool for self-diagnosis or a replacement for medical evaluation, advice, diagnosis or treatment by a healthcare professional. Please speak to your physician if you have questions about your medical condition.

Viewer discretion is advised: Some medical content is graphic and may not be suitable for all viewers.

Attribution Key

for more information see: <http://open.umich.edu/wiki/AttributionPolicy>

Use + Share + Adapt

{ Content the copyright holder, author, or law permits you to use, share and adapt. }



Public Domain – Government: Works that are produced by the U.S. Government. (17 USC § 105)



Public Domain – Expired: Works that are no longer protected due to an expired copyright term.



Public Domain – Self Dedicated: Works that a copyright holder has dedicated to the public domain.



Creative Commons – Zero Waiver



Creative Commons – Attribution License



Creative Commons – Attribution Share Alike License



Creative Commons – Attribution Noncommercial License



Creative Commons – Attribution Noncommercial Share Alike License



GNU – Free Documentation License

Make Your Own Assessment

{ Content Open.Michigan believes can be used, shared, and adapted because it is ineligible for copyright. }



Public Domain – Ineligible: Works that are ineligible for copyright protection in the U.S. (17 USC § 102(b)) *laws in your jurisdiction may differ

{ Content Open.Michigan has used under a Fair Use determination. }



Fair Use: Use of works that is determined to be Fair consistent with the U.S. Copyright Act. (17 USC § 107) *laws in your jurisdiction may differ


Our determination **DOES NOT** mean that all uses of this 3rd-party content are Fair Uses and we **DO NOT** guarantee that your use of the content is Fair.

To use this content you should **do your own independent analysis** to determine whether or not your use will be Fair.


ch216sp12syllabus: Experiment 1

Experiment 1: Synthesis of an Analgesic from Aniline

Introduction: Acetanilide is an analgesic, which was formally known as Antifebrin¹, and is structurally similar to acetaminophen (or Tylenol). However, unlike acetaminophen, acetanilide is toxic. Acetanilide is prepared from aniline using an acetylation reaction. Acetylation is often used to place an acetyl protecting group on primary or secondary amines to reduce their reactivity toward oxidizing agents or electrophiles. Acetamides are usually crystalline solids which can be a help in purification by recrystallization. The melting points can be used for characterization and identification of the corresponding compounds.



acetanilide





acetaminophen (Tylenol)

PD-INEL

1. Cahn, A; Hepp, P. (1886), "Das Antifebrin, ein neues Fiebermittel", *Centralbl. Klin. Med.* **1886**, 7, 561-65

Reaction: Acetylation of aniline and unknown substituted anilines with acetic anhydride



X = unknown substituent group

PD-INEL

Techniques: the [Beilstein test](#), [melting point](#), [recrystallization](#), [gravity](#) & [vacuum filtration](#), [thin layer chromatography](#), and infrared spectroscopy.

Objective: In part 1 you will convert aniline to acetanilide using an acetylation reaction described below. In part 2 you will be given an unknown aromatic amine from table 1-1 below. The amine will be converted into its acetamide analog also using the acetylation procedure described below. Your objective is to compare the two reactions and determine the identity of the unknown amine/acetamide product.

Procedure:

Reagents:

Aniline

HCl (concentrated HCl is 37% w/w)

Acetic anhydride
Sodium acetate

Part 1: Dissolve 500 mg of aniline in 14 mL of water. Note that aniline is immiscible in water and two layers should be observed. Add 0.45 mL of concentrated hydrochloric acid. Measure out 0.6 mL of acetic anhydride and prepare a solution of 530 mg of sodium acetate in 3 mL of water. Add the acetic anhydride to the solution of aniline hydrochloride in water, mix by swirling, and immediately add the sodium acetate solution. The solution becomes white as acetanilide precipitates. Cool the solution in an ice bath and collect the solid acetanilide by vacuum filtration. Recrystallize from 95% ethanol - it may be necessary to use a small amount of water.

Part 2: Use the above procedure for acetylation of the substituted unknown aniline. Exact quantities cannot be calculated since the starting material is an unknown.

Characterization: Both products of part 1 and 2 should be characterized by TLC. Each product should be co-spotted against the starting material for each reaction, respectively, and against each other. Collect the Infrared spectrum of the part 1 product, acetanilide, and the unknown starting material and product from part 2. Characterize the products of part 1 and 2 by m.p and using the Beilstein test for halogens.

For the [Beilstein test](#) first clean a copper wire by holding it briefly in a flame. Touch the wire to a sample of the compound and return it to the flame. A blue-green color in the flame indicates the presence of a halogen. It is a good idea to test the procedure on a compound known to contain a halogen before trying it on your own compound.

Use [Reaxsys](#) or [SDBS](#) to find literature IR spectra of aniline, acetanilide, unknown substituted amine and its corresponding acetamide product for comparison to the spectra that you obtain in lab.

Post-lab Questions (write your answer at the end of your lab notebook pages for this experiment):

1. Describe in your own words what happens to aniline as concentrated hydrochloric acid is added.
2. Why is sodium acetate used? What if NaOH was used instead?
3. Compare the behavior of the unknown substituted amine with aniline in the acetylation reaction. How did the 2 reactions differ? What how might the structural differences in the unknown substituted amine cause it to react differently than aniline, which is un-substituted?
4. How does the Infrared spectrum of your unknown amine compare to that of aniline? How does the spectrum of your unknown acetamide product compare to that of acetanilide?

Table 1-1

Amine	bp (°C)	mp (°C)	Acetamide mp (°C)
aniline	184	--	114
N-methylaniline	196	--	102
<i>ortho</i> -toluidine	200	--	110
<i>meta</i> -toluidine	203	--	65
2-chloroaniline	209	--	87
2-ethylaniline	210	--	111

2,5-dimethylaniline	213	14	139
2,6-dimethylaniline	215	11	177
4-ethylaniline	216	-6	94
2,4-dimethylaniline	217	--	133
3,5-dimethylaniline	220	10	144
2,3-dimethylaniline	221	4	135
3-chloroaniline	230	--	72, 78
<i>meta</i> -anisidine	251	--	81
<i>para</i> -toluidine	200	44	147, 153
3,4-dimethylaniline	224	49	99
2,5-dichloroaniline	251	50	132
<i>para</i> -anisidine	240	58	127
<i>para</i> -bromoaniline	245	66	168
<i>para</i> -chloroaniline	232	72	172, 179